

A Highly Discriminating Taq polymerase Ideal for allele-specific PCR (AS-PCR), SNP detection and CRISPR/Cas gene editing

HiDi® 2x PCR Master Mix Cat no. #9101S, #9101M

Medix
Biochemica

AMP 2022 Annual Meeting & Expo
November 1-5, 2022
Phoenix, Arizona USA

Audrey Crousilles, Laura Vaahtoranta, Ramon Kranaster, Kevin Eboigbodin

Introduction

HiDi® stands for High Discrimination of nucleotides mismatches at the 3'-terminus of PCR primers. The HiDi® 2X PCR master mix is a SYBR Green-compatible ready to use mix that contains a highly selective recombinant DNA polymerase suitable for allele-specific PCR (AS-PCR)¹, primer extension or methylation-specific PCR. HiDi® is similarly effective for checking the quality of genome editing techniques such as CRISPR/Cas or TALEN-based approaches^{2,3}. The enzyme uses a hot-start aptamer-based technology while the optimized master mix chemistry promotes high sensitivity and robust polymerase processivity. Temperatures above 50°-55°C allows the aptamer's secondary structure to melt and will set-free the DNA polymerase. HiDi® efficiently discriminates transversion mutations as well as transition mutations in a more delayed fashion. For probe-based assays, the HiDi® Taq 2X PCR master mix is also available (Cat no. #4200S or #4200M).

Methods

HiDi® 2x PCR Master Mix and three competitors were used to amplify a 100 bp-long fragment of actin gene (Genbank NM_001101.5). A 10-fold dilution series of human genomic DNA was amplified using SYBR Green chemistry over 40 cycles as per Table 1. Matching versus 3'-mismatching reverse primers were designed according to Figure 1. Reaction setup and final primer concentrations were applied according to manufacturer recommendations.

Step	Temperature	Time	Cycles
1	95°C	2 min	1
2	95°C	15 sec	40
3	60°C	30 sec	
4	72°C	30 sec	

Table 1. Thermal cycling conditions. Step 1 of initial denaturation was extended to 10 min for competitor C.

NM_001101.5:85-1212	3' -GTCTTTGCGGATGTCACAGTCACACTTCATGATGGAGT-5'
WT primer	5' -GGATGTCCACGTCCACACTTC-3'
A-mismatch primer	5' -GGATGTCCACGTCCACACTTA-3'
T-mismatch primer	5' -GGATGTCCACGTCCACACTTT-3'
Forward primer	5' -CACTCTTCCAGCCTTCCTC-3'

Figure 1. Primer design. The mismatching reverse primers exhibit either a transversion (A-mismatch primer) or transition (T-mismatch primer) mutation vs. the matching WT primer. All primers have a T_m of ~63°C. WT stands for Wild Type.

Results

Firstly, HiDi® 2x PCR Master Mix outperformed and was the most efficient for discriminating against a transversion mutation (A-mismatch primer) in comparison with the WT primer as shown in panels A and E of Figure 2. The product of the T-mismatch primer is significantly delayed and hence demonstrates good discrimination. Furthermore, the competitors' products yielded lower fluorescence signals and result in poor discrimination of either the A-mismatch or T-mismatch primer as shown in panels B, C, D and E of Figure 2.

Secondly, the efficiency of HiDi® 2x PCR Master Mix after amplification with the matching WT primer stands at 97% and proves to be superior and over the three other products tested for three major suppliers.




Figure 2. Amplification curves and C_q values plots. **A)** HiDi® 2x PCR Master Mix discriminates primers with a 3'-end mismatch. **B)** Competitor A does not discriminate either of the mismatches. **C)** Competitor B performances result in poor discrimination. **D)** Competitor C performances result in poorer discrimination and signals.

E) The C_q values for HiDi® 2x PCR Master Mix and competitors across DNA copies demonstrates the high discriminative power of HiDi® in presence of a transversion or transition mutation.

Note: All three amplification curves presented for WT, A-mismatch and G-mismatch primers correspond to the 50000, 5000 and 500 DNA copy numbers in the samples.

Conclusions

- HiDi® 2x PCR Master Mix displays superior capacity for discriminating against a single mutated nucleotide. It is a product of choice for mutation detection using AS-PCR and for sensitive genome editing techniques.
- HiDi® 2x PCR Master Mix demonstrates superior performance with SYBR Green chemistry compared to three mixes from leading suppliers.
- HiDi® 2x PCR Master Mix is an equally efficient mix for highly sensitive PCR amplification.

References:

1. Miotto, O. et al. Emergence of artemisinin-resistant *Plasmodium falciparum* with kelch13 C580Y mutations on the island of New Guinea. *PLoS Pathog.* 16, 1–21 (2020).
2. Morisaka, H. et al. CRISPR-Cas3 induces broad and unidirectional genome editing in human cells. *Nat. Commun.* 10, (2019).
3. Sakurai, T. et al. Bindel-PCR: a novel and convenient method for identifying CRISPR/Cas9-induced biallelic mutants through modified PCR using *Thermus aquaticus* DNA polymerase. *Sci. Rep.* 9, 1–14 (2019).

Medix Biochemica Group

Kloivpellontie 3, FI-02180 Espoo, Finland

medix@medixbiochemica.com • www.medixbiochemica.com