

A close-up photograph of a laboratory procedure. A white pipette with a blue tip is being used to dispense a clear liquid into one of the wells of a green microplate. The pipette is held by a gloved hand. The background is a blurred laboratory environment with blue chairs and a white bench. The text "Novel Engineered Polymerases For The Detection of Infectious Disease" is overlaid in the center of the image.

Novel Engineered Polymerases For The Detection of Infectious Disease

Medix Biochemicals

Before we get started...



Please use the Q&A button at the bottom of your screen to submit any questions.

There will be a 10-minute Q&A session at the end of this webinar. If we do not have time to answer your questions, we will reach out by email following the webinar.



This webinar is being recorded. The recording as well as on-demand link will be released following the webinar.


Medix Biochemica

We're Wherever You Are

- A truly global company with a proudly Finnish background
- Our products are used in over 70 countries
- Offices and manufacturing facilities all around the world

Employee Expertise

- Decades of dedication to developing essential raw materials for diagnostic applications
- Approximately 270 employees, with 1 in 5 people work in R&D



Engineered Polymerases for Molecular Diagnostics and more...

- Unique Polymerases for DNA/RNA
- ISO:13485 Conformity
- Assay Development Services
- Lyophilization Services
- Lyo-Ready PCR Products



Today's Speaker:

**Dr. Maja
Studencka-Turski**
Scientific Lead R&D



01


Engineered DNA Polymerases in Biotechnology

Medix Biochemica

DNA Polymerases Play Central Roles In Modern Molecular Biology And Biotechnology




DNA Cloning



DNA Sequencing




Molecular Diagnostics



Polymerase Chain Reaction (PCR)



Synthetic Biology




SNP Detection

DNA Polymerases – “Nature’s Molecular Machines”

Each thermostable DNA polymerase has its own set of unique characteristics:

- Thermostability
- Extension rate
- Fidelity
- Processivity
- Specificity
- Ability to bypass damage
- Nuclease activity
- Strand displacement activity





“


“Natural DNA Polymerases Often
Do Not Have The ‘Performance
Specifications’ Needed For
Transformative Technologies”

(Zahra Ouaray, *Et Al.* J Biol Chem. 2020 Dec 11; 295(50): 17046–17059.)

What does Engineering Polymerases Allow Us To Do?


Each thermostable DNA polymerase has its own set of unique characteristics:

- Thermostability
 - Extension rate
 - Fidelity
 - Processivity
 - Specificity
 - Ability to bypass damage
 - Nuclease activity
 - Strand displacement activity
- Engineering Polymerases allows us to adjust any of these characteristics to solve a wide range of problems**



How Do We Genetically Engineer DNA Polymerases?

Tailoring DNA Polymerases By **Random Mutagenesis**



Rational Design Approach

- Sufficient Information About The Enzyme Necessary
 - e.g. Site Directed Mutagenesis
- Example: KlenTaq Polymerase → Mutation Of Amino Acids Involved In Substrate Binding



02

Infectious Disease Applications

Medix Biochemica

Engineered Polymerases And Their Role In Infectious Disease

Engineering polymerases to have unique characteristics could significantly advance the field of infectious disease testing.

Improvements and innovations could include:



Increased Sensitivity and Specificity

- More sensitive detection of pathogens, (i.e. lower limit of detection)
- Enhanced specificity to reduce likelihood of false positives



Faster Reaction Times

- More efficient polymerases to speed up the polymerase chain reaction (PCR) process.
 - Crucial in infectious disease diagnostics where rapid results can lead to quicker treatment decisions and containment measures
 - Reliable results with crude samples without extraction step



Broader Range of Detectable Pathogens

- Engineered polymerases could recognize a wider range of sequences or to function under varying conditions
- Tests could be developed to detect a broader spectrum of pathogens, including emerging or mutating infectious agents.

Engineered Polymerases And Their Role In Infectious Disease

Improved Stability and Versatility

- Engineered polymerases could be more stable at room temperature, or function effectively in a wider range of conditions
 - Making diagnostic tests **more versatile and easier to distribute**.
 - Particularly important in resource-limited settings or **in-field diagnostics** where laboratory facilities are not available.
 - Key in enabling **Point-of-Care testing**
- As well as enabling the **Integration with Advanced technologies** such as NGS and CRISPR-based technologies



Our Engineered Polymerases



PlexTaq® 5x qPCR
Multiplex Master Mix



RNA Applications
(RT-PCR)



Isothermal Amplification

PlexTaq® 5x qPCR Multiplex Master Mix



Medix Biochemica

PlexTaq[®] 5x qPCR Multiplex Master Mix

All-In-One - A Unique Engineered Solution

Robust - Up to 30 target multiplexing in real-time

5X Concentration - Maximizing volume available for primers/probes



Fast Time-To-Results - Reliable results with crude samples (High Inhibitor Tolerance)

Highly Specific - No false amplification.

Lyo-ready - Enables RT storage and shipping once dried.

DirectPCR From Blood

– PlexTaq[®] 5x qPCR Multiplex Master Mix



Traditional Approach

Blood Samples

DNA Extraction

qPCR

Direct Approach

Blood Samples

Direct to qPCR

qPCR



No extraction reagents or equipment required



Less error-prone




Higher throughput

Real-time DirectPCR from blood

– PlexTaq[®] 5x qPCR Multiplex Master Mix




K3 EDTA blood




FAM channel: Human coagulation risk factor 5- SNP gene target

HEX channel: spiked internal control target (10^4 c/rxn)

Cy5 channel: ACTB1 target (human gene present in human blood sample)



Agarose gel after DirectqPCR from blood – PlexTaq[®] 5x qPCR Multiplex Master Mix




Equally good and specific target amplification regardless of blood concentration.

Multiplex Microarray

– PlexTaq[®] 5x qPCR Multiplex Master Mix

Mycobacterium tuberculosis detection by microarray

Spotting scheme of the human β -actin, KRAS and RRDR multiplex




ISTITUTO DI SCIENZE E TECNOLOGIE CHIMICHE GIULIO NATTA

SCITEC National Research
Council – CNR, Italy

Multiplex Microarray – PlexTaq[®] 5x qPCR Multiplex

Fluorescence intensities for the three multiplex targets



PlexTaq® 5x qPCR Multiplex Master Mix

Conclusion

- PlexTaq® 5x qPCR Multiplex Master Mix is a **5X concentration** engineered polymerase, **ideal for multiplexing**.
 - Widely used “workhorse” mix
 - Up to **30 target** multiplexing
 - **Lyo-ready** formulation for lyophilization
 - Highest quality standard - **ISO 13485:2016** conformity
- Direct real-time qPCRs on crude samples – **No Nucleic Acid Extraction**
- Real time multiplex detection qPCRs – **Food Pathogen Testing**

RNA Applications (RT-PCR)



Medix Biochemica

Volcano® RT-PCR Master Mix

Engineered Taq DNA polymerases with reverse transcriptase activity

Unique - Simultaneous reverse transcription and amplification

Ready-To-Use - Simple and sensitive


Fast Time-To-Results - Reliable results with crude samples

Hot-Start - Aptamer formulation prevents unspecific amplification at lower temperatures

Thermostable – Active at temperatures where complex samples (e.g. secondary RNA structures) are broken down

Volcano[®] RT-PCR Master Mix – Thermostability

Unique - Truly **HOT** Reverse Transcription



Performing reverse transcription at a higher temperature can be beneficial for assays dealing with complex samples

Blatter, N., Bergen, K., Nolte, O., Welte, W., Diederichs, K., Mayer, J., Wieland, M. and Marx, A. (2013), Structure and Function of an RNA-Reading Thermostable DNA Polymerase. *Angew. Chem. Int. Ed.*, 52: 11935–11939. doi: 10.1002/anie.201306655


R. Kranaster, M. Drum, N. Engel, M. Weidmann, F. T. Hufert and A. Marx, (2010), One-step RNA pathogen detection with reverse transcriptase activity of a mutated thermostable *Thermus aquaticus* DNA polymerase, *Biotechnol. J.*, 5(2), 224-31.

Kranaster R, Zeller J, Kühn B, Marx A, (2016) Neues Enzym mit Reverse Transkriptase- und DNA-Polymerase-Funktion, *BioSpektrum*, 2, 164-165.


Volcano[®] RT-PCR Master Mix – 0-Step RT-qPCR

Fast start function - 0-step RT-PCR

Traditional 1-Step Approach



0-Step RT-PCR (Volcano)



Simplified RT-PCR without isothermal RT-steps



Compatible with dsDNA-sensitive dyes (e.g. SYBR)



Additional polymerases can be added

Medix Biochemica

Volcano® RT-PCR Master Mix – DirectPCR

Wastewater screening of viruses

1. SARS-CoV-2
2. HIV
3. Monkeypox

“(…) the entire process was accelerated by a new analytical technique called V2G (for “volcano second generation”) quantitative polymerase chain reaction (or qPCR) that can pick up COVID-19 in sewage faster than earlier methods.”
(after University of Miami))

Qualitative and quantitative detection of viral targets:

1. HIV-1
2. Zika virus



Contents lists available at ScienceDirect
Science of the Total Environment
journal homepage: www.elsevier.com/locate/scitotenv

ELSEVIER

Lessons learned from SARS-CoV-2 measurements in wastewater
Mark E. Sharkey^a, Naresh Kumar^b, Alejandro M.A. Mantero^b, Kristina M. Babler^c, Melinda M. Boone^d, Yoslayma Cardente^d, Elena M. Cortizas^d, George S. Grills^d, James Herrin^e, Jenny M. Kemper^d, Richard J. Christopherson^d, Natasha J. Taubert^d, John J. Tamplin^d, Benjamin J. Gold^d

First detection of the Monkeypox virus using wastewater-based surveillance in Miami-Dade County

Mark E Sharkey, Kristina M Babler, Ayaaz Amirali, George S Grills, and 10 more

This is a preprint; it has not been peer reviewed by a journal.

RESEARCH ARTICLE
Detection of SARS-CoV-2 from raw patient samples by coupled high temperature reverse transcription and amplification

Johannes W. P. Kuiper^{1,2}, Timo Baade^{1,2}, Marcel Kremer³, Ramon Kranaster⁴, Linda Irmisch⁵, Marcus Schuchmann⁵, Johannes Zander³, Andreas Marx^{2,4,6}, Christof R. Hauck^{1,2*}

MedixMDx qRT-PCR Lyo ready Mix

Sensitive and Versatile RT-PCR mix ready for lyophilization

Lyo-Ready Mix

- Specifically designed for lyophilization. No need for excipients.

Robust And Sensitive

- Optimized for rapid detection and quantification of various RNA templates, including mRNA, viral RNA, and total RNA.
- Effective for detecting multiplexed viral targets directly from human patient samples.


Versatile

- Antibody-regulated hot-start Taq polymerase and reverse transcriptase, along with target-specific primers and probes, and RTase Lyo.

One-Step Probe Mix

- Universal one-step probe mix.
- Directly lyophilized without additional components.

MedixMDx qRT-PCR Lyo Ready Mix



RNA Applications (RT-PCR)

Conclusion

Volcano3G

- **Simultaneous** reverse transcription and amplification
- **Simple** and **sensitive**
- **Reliable** results with crude samples
- Aptamer-mediated **hot-start** - Active at temperatures where **complex sample structures** are broken down
- Highest quality standard - **ISO 13485:2016** conformity
- **Rapid detection and identification of RNA & DNA targets** (e.g., SARS-CoV-2, RSV and Flu A/B screening)
- **Lyo-ready** formulation for lyophilization available upon request

MedixMDx qRT-PCR Lyo Ready Mix

- **Sensitive**
- **Robust**
- **Lyo-Ready**
- **Multiplex Capabilities**

Isothermal Amplification



Medix Biochemica

Isothermal Amplification

Range of Engineered Bst Polymerases

MedixMDx Fast Bst Product Line

- Perfectly suited for Isothermal Amplifications
- Options for both LAMP and RT-LAMP
- Fast time-to-results
- Available as Master Mixes or separated components

Isotherm3G Product Line

- Single-enzyme LAMP/RT-LAMP
- POC applications

Perfectly Suited For Loop-Mediated Isothermal Amplification (LAMP)

Loop-mediated isothermal amplification (LAMP) is a single-tube technique for the amplification of DNA and a **low-cost** alternative to detect certain diseases.



Reverse transcription loop-mediated isothermal amplification (RT-LAMP) combines LAMP with a reverse transcription step to allow the detection of RNA.

Advantages of LAMP in Molecular Diagnostics



- ✓ No Dedicated Equipment Needed
- ✓ Fast Reaction (~ 15 min)
- ✓ Easy Visualization
- ✓ Simplicity of Operation
- ✓ Ideal for Point-of-Care or At-Home Use



MedixMDX – Fast Bst Family




Faster Time-To-Results without sacrificing performance




High Inhibitor Tolerance

Warm-Start BST + UDG


Target: *Salmonella enterica*



MedixMDx Fast Bst



MedixMDx Warm-Start Bst
+ UDG



Competitor (Warm-Start)
+ UDG

Collaboration with:





ISTITUTO
ITALIANO DI
TECNOLOGIA

Medix Biochemica

Isotherm3G with Improved RT-Activity



- Fast and specific detection
- No addition of a reverse transcriptase is needed!





Reliable detection of SARS-CoV-2 N gene

Isotherm3G vs Competitor

Isotherm3G



Competitor



Salmonella

Similar results when also tested on:
Staphylococcus, and Pseudomonas

Legionella

Collaboration with:



ISTITUTO
ITALIANO DI
TECNOLOGIA

Isothermal Amplification

Conclusion

MedixMDx Fast Bst & Isotherm3G Product Lines

- Ideal for **Point-of-Care** Applications
- Products for both **LAMP** and **RT-LAMP**
- **Fast Time-To-Results**
- **Simultaneous reverse transcription and amplification** with Isotherm3G
- Available as Master Mixes or separated components
- Highest quality standard - **ISO 13485:2016** conformity

Q&A Session



Dr. Maja Studencka-Turski



Dr. Arielle Bryan

A close-up photograph of two hands shaking in a firm grip. The scene is lit with vibrant blue and red neon lights, creating a modern, high-tech atmosphere. A thin white circle is superimposed over the handshake, framing the central text.

Thank you

Medix Biochemica

Don't Miss Out on Our Next Webinars

Subscribe to our events:



Follow Medix Biochemica on social media:

