


A person wearing a blue lab coat and gloves is holding a test tube containing a green liquid. The test tube is positioned over a multi-well plate, which also contains several wells with green liquid. The background is a blurred laboratory setting. The entire image has a blue tint.

Improving Point-of-Need Detection with LAMP:

FROM SNP ANALYSIS TO SENSITIVE PATHOGEN DETECTION

Before we get started

 Please use the Q&A button at the bottom of your screen to submit any questions.

There will be a Q&A session at the end of this webinar. If we do not have time to answer your questions, we will reach out by email following the webinar.

 This webinar is being recorded. The recording as well as on-demand link will be released following the webinar.

Independent, International, and Industry-Leading Raw Materials Supplier

- Provider of high quality antibodies, antigens, proteins, enzymes
- Expertise in immunoassays, clinical chemistry
- Engineered Polymerases for Molecular Diagnostics
 - Unique Polymerases for DNA/RNA
 - ISO:13485 Conformity
 - Assay Development Services
 - Lyophilization Services
 - Lyo-Ready PCR Products



Featured Speakers and Panelists



Dr. Paola Cecere

Postdoctoral Researcher
Istituto Italiano Di Tecnologia



Dr. Pier Paolo Pompa



Research Director
Istituto Italiano Di Tecnologia



Dr. Giuseppina Sannino

Global Business Development Manager, MDx
Medix Biochemica





Improving Point-of-Need Detection with LAMP: from SNP Analysis to Sensitive Pathogen Detection

Dr. Paola Cecere
Postdoctoral Researcher – Istituto Italiano Di Tecnologia

29th October 2024



Dr. Pier Paolo Pompa

Principal Investigator

Nanodiagnostic platform



Integrated, smart, low-cost, and rapid assays/sensors for on field and point-of-care diagnostics.

Specific areas of applications span from clinical diagnostics (cancer, genetic, and infectious diseases) to food safety and traceability, and environmental control (pathogens, pollutants and contaminants).

This research line is based on hybrid detection strategies:

- green and controlled synthesis of nanomaterials
- surface (bio)chemistry
- plasmon non-linear response
- nanozymes
- biotechnology and molecular biology



**LAMP-based
GENETIC TESTING**

LAMP: Loop-mediated isothermal AMplification

© 2000 Oxford University Press

Nucleic Acids Research, 2000, Vol. 28, No. 12

e63


Loop-mediated isothermal amplification of DNA

Tsugunori Notomi^{1,3,*}, Hiroto Okayama², Harumi Masubuchi¹, Toshihiro Yonekawa¹,
Keiko Watanabe¹, Nobuyuki Amino³ and Tetsu Hase¹

¹Eiken Chemical Co. Ltd, 1381-3 Shimoishigami, Ohtawara, Tochigi 324-0036, Japan, ²Department of Biochemistry and Molecular Biology, The University of Tokyo, Graduate School of Medicine, Bunkyo-ku, Tokyo 113-0033, Japan and ³Department of Laboratory Medicine, Osaka University Medical School, 2-2 Yamadaoka, Suita, Osaka 565-0871, Japan

- Outstanding reaction for DNA and RNA amplification, with broad and multifaceted capabilities
- Robustness and accessibility

Mastering
design to
maximize the
elevated LAMP
performances



Pirelli's iconic slogan

Our Solution:



Medix Biochemica

“Lyo-LAMP”

Unleashing LAMP potential with our expertise!

Section 1

1. How does LAMP work?
2. What are LAMP's strengths?
3. Which challenges do we face with LAMP?
4. Where can LAMP be carried out?

Section 2


Our answer: "Lyo-LAMP"

1. Lyo-LAMP identity card
2. Deep dive into Lyo-LAMP
3. Applications and perspectives

The answer is in the name


PCR

(Polymerase Chain Reaction, Kary Mullis, 1983)




thermal cycles

LAMP



Loop-mediated



Loop: ssDNA

Stem: dsDNA

isothermal amplification




single reaction temperature
(60-65 °C), no denaturation
step

Exploring the LAMP mechanism – things to keep in mind

LAMP primers
set of 4-6 specially designed primers, specific for 6-8 target regions

- FIP = forward inner primer
BIP = backward inner primer
 - F3 = forward outer primer
B3 = backward outer primer
 - FL = forward loop primer
BL = backward loop primer
- Guide for the formation of the loop-stem structures
- Work only in stage 1
- Accelerating primers, work only in stage 2



Exploring the LAMP mechanism – things to keep in mind

Alternative mechanisms to the denaturation stage in LAMP


1. dynamic equilibrium of the dsDNA target at 60-65 °C
2. DNA polymerase with high strand-displacement activity
3. two priming sites for each stem-loop structure:



Loop: ssDNA

Stem: dsDNA


a. Primer annealing



b. Self-priming


STAGE 1: production of the starting material

STAGE 2: auto-cycling strand-displacement amplification/elongation step




STAGE 1: production of the starting material

4. F3 annealing




5. F3 elongation



6. polymerase
strand displacement
activity

STAGE 1: production of the starting material



STAGE 1: production of the starting material




10. B3 annealing

11. elongation

12. BIP-elongated strand
ssDNA displacement

13. loop folding at both ends



DUMBELL-LIKE STRUCTURE


DUMBELL-LIKE STRUCTURE



artificial ssDNA target

STAGE 2: auto-cycling strand-displacement amplification/elongation step


STAGE 2: auto-cycling strand-displacement amplification/elongation step



MULTIPLE
AMPLIFICATION STARTING POINTS
on the dumbbell-like structures:

1. annealing inner primer/loop
2. self-priming (3' stem end)
3. annealing loop primer/loop

STAGE 2: auto-cycling strand-displacement amplification/elongation step



- Multiple amplification sites
 - Loop unfolding-refolding
 - Strand-displacement activity
-
- Alternatively inverted repetitions of the target
 - Cauliflower-like structures


**Exponential
Auto-cycling
Amplification**

Agarose gel



What are LAMP's strengths?

	PCR the precision powerhouse for lab-based quantification	LAMP the speedy and hassle-free tool for on-the-spot diagnostics	
Temperature	Thermal cycles	Isothermal	
Equipment	Expensive thermal cycler	Minimal instrumentation	→ Cost-effective
Robustness	Sensitive to inhibitor	Highly tolerant to inhibitors	→ Simplified extraction protocols
Amplification efficiency	10 ⁶	10 ⁹ – 10 ¹⁰	→ Greater yield



What are LAMP's strengths?

	PCR the precision powerhouse for lab-based quantification	LAMP the speedy and hassle-free tool for on-the-spot diagnostics	
Temperature	Thermal cycles	Isothermal	
Equipment	Expensive thermal cycler	Minimal instrumentation	→ Cost-effective
Robustness	Sensitive to inhibitor	Highly tolerant to inhibitors	→ Simplified extraction protocols
Amplification efficiency	10^6	$10^9 - 10^{10}$	→ Greater yield
		High-performing LAMP	
SPEED	> 1 hour	7 - 40 mins	→ Faster time-to-results
SENSITIVITY	10-100 copies/reaction	1-3 copies/reaction	→ Greater sensitivity
SPECIFICITY	2 primers x 2 target regions	4-6 primers x 6-8 target regions	→ Greater specificity

Which challenges do we face with high-performing LAMP?

1. Complex primer design

Overall assay performance:

- sensitivity (false negative %)
- specificity (false positive %)
- reaction rate
- reliability/reproducibility

Detection of challenging DNA markers (as SNPs)

- Chemical and thermodynamic requirements
- Specific detection
- Undetectable background
- Avoid non-specific binding and primer dimers

**Need for
extensive
expertise in
primer design**

2. Mix composition

- Overall assay performance
- Efficient warm-start polymerase
- Lyo-compatible

3. Risk for carryover contaminations

The assay must include a system against cross-contaminations

4. Extraction protocol

Simple, fast and LAMP-compatible

A **high-performing LAMP** can be particularly suitable for

Genetic testing

Point-of-care

Rapid testing at the point of need (centralized or decentralized settings)

- Fast turnaround time
- User-friendly
- Affordability
- Portability

High-throughput

Continuous and large-scale screenings

- Scalability
- High sample capacity
- Automation
- Reproducibility

Main requirements for POC and HT testing design:

- **Chemistry:** fast and reliable, no cross-contaminations, inhibitor tolerance
- **Workflow:** short and easy
- **Supply Chain:** low-cost components with high performance, scalability, room storage
- **Instrumentation:** wide equipment availability, cost-effectiveness



Expertise in LAMP technology



Medix Biochemical


Expertise in enzyme engineering
and lyophilization

Our Answer:
“Lyo-LAMP”

High-performing LAMP
in Lyophilized format
for molecular diagnostic solutions



Smart



User-Friendly




Business-Friendly



Smart

- All-inclusive solutions → DNA extraction & amplification
- High-performing LAMP →
 - Customized primer set
 - High sensitivity and specificity
 - Fast results
- Lyophilized format → RT-shipping and storage, with long-term stability



User-Friendly

- Highly tolerable test
- Easy-to-use
- Minimal instrumentation required



Business-Friendly

- Cost-effective
- Alternative solutions to the limits of traditional methods
- Scalable
- Custom on-demand

“Lyo-LAMP” is suitable for the detection of different **DNA markers**, such as:


- Species-specific genes



- SNP-genotyping
(Single nucleotide polymorphism)



- InDel-genotyping
(insertion-deletion)



Targets:

- Human genome
- Farm-animal and pet genome
- Microorganism genome (bacteria, parasites, virus)
- Food

Applications

- Infectious diseases
- Genetic predisposition and susceptibility
- Pharmacogenetics and nutrigenetics
- Food safety and traceability
- Environmental monitoring (water, biosolids, waste-water)
- Bio-terrorism

1

Sampling

Different (target-depending)
sample types,
non-invasive procedures

Human and veterinarian samples



saliva – buccal/nasal swab - fingerprick

Food and beverages samples



1

Sampling

2

**Streamlined protocol
for DNA extraction**

Thermo-chemical
(sample-depending)
DNA extraction

- Ready-to-use extraction solution with RT storage
- Simple heater required
- A few simple operating steps



No DNA purification!



Total processing time: 5-10 mins!

1

Sampling

2

Streamlined protocol
for DNA extraction

3

LAMP amplification
step

Total reaction
time
from 7 to 40 mins!

- High-performing LAMP with specific primers
- LyoBeads technology
- Fluorescence-based detection (real-time)
- Open test (no test-associated instrument)
- Simple instrumentation required


The test design prevents **LAMP carry-over cross-contaminations**:

- ✓ Using a one-pot closed-tube format for amplification and detection;
- ✓ Integrating an enzymatic system for contamination degradation.

“save-time” test
for fast screening




Deep dive into “Lyo-LAMP” – Customized lyo-bead technology



- Pre-dispensed LAMP mix in freeze-dried spheres (1 bead = 1 reaction)
- Specific customized primer set design


Simple handling




- Ready-to-use
- Customized for our tests (optimized bead composition)
- Long-term stability, room temperature shipping and storage
- Scale-up production

Salmonella enterica

Competitor in solution




Lyo-LAMP



1

Real-Time SNP and InDel-Genotyping




SNPs (single nucleotide polymorphism)

genetic point mutations, that recur with a frequency of at least 1% in a species population

InDels (insertion-deletion polymorphism)

indels can be small, involving just a few base pairs, or larger, spanning several base pairs



Ideal marker

SNPs are involved in development and progression of different diseases, in the genetic susceptibility and predisposition to specific conditions, in pharmacogenetics, in varietal discrimination, etc.

Challenging marker

- Complex and expensive procedures
- High background signal
- High-quality DNA samples

Traditional methods


e.g. microarrays, PCR-RFLP, TaqMan assays, Next-generation sequencing...




POC testing



Double-tube test		TARGET GENOTYPE		
		SNP-mutated	SNP-wild type	SNP-heterozygous
 tube 1 mut LAMP-probe		Perfect match ✓	single-base mismatch ✗	Perfect match ✓
 tube 2 wt LAMP-probe		single-base mismatch ✗	Perfect match ✓	Perfect match ✓



Real-Time SNP-Genotyping with “Lyo-LAMP”: The Detection Strategy



3 simple steps!


ALL-INCLUSIVE TEST

Components:


- Ready-to-use extraction solution
- mut lyo-bead tube
- wt lyo-bead tube

Test results:

- homozygous mutation
- heterozygous
- no mutation



Collect the buccal swab




95 °C, 5 min

Place and stir the swab into the extraction solution

Place the tube into a heater set at 95 °C, for 5 minutes

Purification-free

mut wt



Resuspend the wt and mut lyobeads and add the extract into each tube

Place the 2 reaction tubes into the instrument

The test is completed within 30-35 minutes

1. Sample Collection

20 seconds


2. Sample Prep

5 minutes

3. Test Execution

35 minutes

Target	Marker type	Sample	Application
<i>MCM6</i> gene	13910 C>T (SNP)	buccal swab	Lactose intolerance
<i>MTHFR</i> gene	677 C>T (SNP)	buccal swab	Hyperhomocysteinemia and folate level
<i>CYP2C19*2</i>	681 G>A (SNP)	buccal swab	Antiplatelet drug (clopidogrel) metabolism
<i>CYP2C19*3</i>	636 G>A (SNP)	buccal swab	Antiplatelet drug (clopidogrel) metabolism
<i>ACE</i> gene	rs4340 in287 (InDel)	buccal swab	Salt sensitivity
<i>MDR1</i> gene (<i>dog</i>)	nt230 del4 (InDel)	buccal swab	P-gp dependent drugs toxicity



Genetic Predisposition to Lactose Intolerance


Opportunity for customizable targets!

2


Food Fraud and Traceability

Food Fraud and Traceability with “Lyo-LAMP”: Test Panel

Target	DNA marker	Matrix	Sensitivity	Application
Triticum aestivum (common wheat) vs Triricum durum (durum wheat)	species-specific gene	wheat	0.3% adulteration	species discrimination
Durum wheat «aureo» variety vs durum wheat «non aureo» varieties	SNP on Chr7A	wheat	10% adulteration	varietal discrimination



Species Discrimination



Varietal Discrimination




Routine Quality Controls

3

Microorganism Detection


Microorganism detection with “Lyo-LAMP”: Test Panel

Target	Matrix/sample	Application
<i>Escherichia coli</i>	tap water milk	water contamination bovine mastitis
<i>Legionella pneumophila</i>	tap water	water contamination
<i>Aeromonas hydrophila</i>	tap water	water contamination
<i>Salmonella enterica</i>	biosolids	biosolids contamination
<i>Enterococcus faecium</i>	tap/waste water	water contamination
<i>Pseudomonas aeruginosa</i>	milk	bovine mastitis
<i>Streptococcus pneumoniae</i>	nasal swab	sinusitis
<i>Staphylococcus aureus</i> (MRSA vs. MSSA)	nasal swab	sinusitis
<i>Clostridium tyrobutyricum</i> (spores)	milk	cheese spoilage
<i>Mycobacterium tuberculosis</i>	saliva	tuberculosis
<i>Plasmodium spp.</i>	finger prick	malaria




Reaction time: 10-20 mins
Sensitivity: 2-10 copies/reaction

Mycobacterium tuberculosis



Salmonella enterica



Opportunity for customizable targets!

Traditional Methods Include:

SMEAR microscopy

Culture

PCR based methods

Antigen detection


Phage based assays

Tuberculin skin testing



Gold-standard POC competitor (example)

- PCR based
- ~77 minutes TTR
- Device-associated test



ALL-INCLUSIVE TEST

Components:

- ready-to-use extraction solution;
- reaction tube with the LAMP-lyobead


Test results:

- Positive or Negative

3 simple steps!



Collect the saliva sample




95 °C, 5 min

Add saliva into the extraction solution;

Place the tube into a heater set at 95 °C, for 5 minutes.

Purification-free



Resuspend the lyobead and add the extract

Place the reaction tube into the instrument

The test is completed within 10-20 minutes

1. Sample Collection

20 seconds

2. Sample Prep


5 minutes

3. Test Execution

10-20 minutes

Lyo-LAMP

High-performing LAMP in Lyophilized format



C
u
s
t
o
m
i
z
a
b
l
e

SAMPLE matrix and volume

Technology FLEXIBILITY

- Challenging targets (SNP)
- Ultrasensitivity detection (e.g., microorganisms)

Application VERSATILITY

- Infectious diseases
- Genetic predisposition and susceptibility
- Pharmacogenetics and nutrigenetics
- Food safety and traceability
- Environmental monitoring (water, biosolids, waste-water)
- Bio-terrorism

If you are interested in collaborating, we would be excited to hear from you!

Amplify Diagnostics With Lyo-LAMP, your partner in Point-of-Care Excellence!



Dr. Paola Cecere
Postdoctoral Researcher –
Istituto Italiano Di Tecnologia

Dr. Pier Paolo Pompa
Research Director –
Istituto Italiano Di Tecnologia



Dr. Giuseppina Sannino
Global Business Development Manager –
Medix Biochemica

Thank You – Any Questions?




Thank you

Medix Biochemica

Did you enjoy the webinar?

Share your feedback:



Medix Biochemica

**Subscribe to stay informed on all
Medix Biochemica webinars:**



Medix Biochemica